

Facilitated genotoxicity screenings in expanded human upcyte® keratinocytes

Kathleen Jarnouen², Astrid Nörenberg¹, Timo Johannssen¹, Torge Evenburg¹, Meryl Roudaut², Julian Bursztyka², Pierre-Jean Ferron²

¹ upcyte technologies GmbH, Osterfeldstr. 12-14, 22529 Hamburg, Germany

² HCS Pharma, 6 Rue Pierre Joseph Colin, 35000 Rennes, France


INTRODUCTION

Genotoxicity assessment of pharmaceuticals or active ingredients used in cosmetology is mandatory for registration. In addition to the standard *in vitro*/*in vivo* test battery, other assays can be of interest because of their high throughput in the drug discovery stage. Among them, the γ H2AX *in vitro* assay was suggested for detecting genotoxic properties of chemicals. The phosphorylation of γ H2AX has been demonstrated to be a sensitive marker for DNA double strand breaks (DSB). Likewise, p53-binding protein 1 (53BP1) is an important regulator of the cellular response, promoting the end-joining of distal DNA ends.

Keratinocytes are routinely used for skin toxicity studies. However, large-scale cultures of primary keratinocytes demanded by high throughput approaches face obstacles such as a limited lifetime and *in vitro* differentiation, often compromising reproducibility due to inhomogeneous cultures. Here, we describe the expansion of primary keratinocytes after lentiviral transduction (upcyte® keratinocytes) for large-scale production and their subsequent use in γ H2AX/53BP1 genotoxicity screenings.

RESULTS


the upcyte® technology



transduction of primary cells with a defined cocktail of lentiviral vectors

Previously, a novel technique was reported which enabled the generation of non-transformed proliferating liver cells from primary human hepatocytes with maintained adult phenotype, so-called "upcyte® cells". upcyte® cells start to grow from primary cells after transduction with a defined cocktail of lentiviral vectors carrying proliferation inducing genes. upcyte® cells have the ability to proliferate for additional cell doublings, depending on the cell type, without losing functional and phenotypic characteristics of mature cells.

morphology & surface markers of upcyte® keratinocytes



analysis of upcyte® keratinocytes


The morphology of upcyte® keratinocytes is similar to the morphology of primary keratinocytes. upcyte® keratinocytes are analyzed by fluorescence microscopy for the expression of the **cytokeratins 5/14 (red, middle panel), 1/10 (red, right panel)** and counter-stained against **DAPI (blue)**. They are characterized by high expression of the basal marker CK5/14 while lacking expression of the suprabasal markers CK1/10.

suitability of upcyte® keratinocytes for genotoxicity

Mean γ H2AX intranuclei intensity after 24 h

Mean γ H2AX intranuclei spot count after 24 h

Mean 53BP1 intranuclei spot count after 24 h



To further evaluate suitability of upcyte® keratinocytes for genotoxicity studies, upcyte® keratinocytes and neonatal human epidermal keratinocytes (NHEK) were challenged with 16 model compounds including **non-genotoxic controls** (e.g. ampicillin trihydrate), **clastogens** (e.g. mitomycin C (MMC)) and **aneugens** (e.g. paclitaxel).

Intranuclear γ H2AX and 53BP1 spots were analyzed using an automated imaging approach. When looking at the γ H2AX spot count, upcyte® keratinocytes were comparable to NHEK and allowed accurate identification of clastogenic compounds. In contrast, no increase of γ H2AX was observed with non-genotoxic substances.

Interestingly, intranuclei 53BP1 spot counts revealed a base level 2 to 4.8 fold higher in upcyte® keratinocytes when compared to NHEK. This may explain that the fold induction of intranuclear 53BP1 spots after exposure to genotoxic compounds tended to be smaller in upcyte® keratinocytes when compared to NHEK.

SUMMARY & CONCLUSION

In conclusion, fold change of intranuclei γ H2AX spot count in upcyte keratinocytes follow the same trends as in NHEK. On the other hand, intranuclei 53BP1 spot count shows a base level 2 to 4.8 fold higher in upcyte keratinocytes when compared to the batch of NHEK used (DMSO 0.5%). Consequently, upcyte® keratinocytes were suitable for genotoxicity studies using the γ H2AX assay for DSB assessment. Combined with virtually unlimited cell access, they represent an attractive platform for automated screening approaches such as High Content Imaging.