

# Generation of Expanded Primary Cells for Cell-Based Toxicity and Metabolism Screenings

Nikolett Nagy, Torge Evenburg, Sandra Rohmoser, Astrid Noerenberg, Timo Johannssen

upcyte technologies GmbH, Osterfeldstr. 12-14, 22529 Hamburg, Germany


## INTRODUCTION

Cell-based assays are a valuable tool to predict *in vivo* effects of drug candidates during early steps of development. Cell-based assays are performed using either cell lines or primary cells. Most cell lines are easy to handle and offer the advantage of infinite proliferation, allowing the generation of large cell banks and a facilitated use in screenings or long-term experiments. However, due to their transformed phenotype, many cell lines often exhibit a reduced physiological relevance. In contrast, primary cells are more representative of the *in vivo* state when compared to cell lines. However, their use *in vitro* is hampered by limited tissue availability, scarce cell yields and a restriction or even lack of proliferation. Taken together, these factors may significantly compromise the scope, length and reproducibility of experiments and often circumvent their use for extended cell-based screenings.

Here, we describe the controlled expansion of human primary cells by lentiviral transduction with proliferation-inducing genes, enabling production volumes of up to 2500 vials containing  $5 \cdot 10^6$  cells each. As a proof of principle, primary cells from several relevant target tissues (liver, skin, kidney, lung) were transduced, subsequently demonstrating successful expansion to large master and working cell banks.

## RESULTS

### generation of upcyte® cells




#### transduction of primary cells with a defined cocktail of lentiviral vectors

We first generated a library of lentiviral vectors carrying proliferation-inducing genes, allowing primary cells to bypass senescence. Different primary cells such as hepatocytes, liver sinusoidal endothelial cells, keratinocytes, proximal tubular epithelial cells and bronchial epithelial cells were transduced. Resulting upcyte® cells gained the ability to proliferate for up to 40 additional population doublings (PDs) without losing functional and phenotypic characteristics of mature cells. All cells exhibited expected morphology patterns and were restricted by the presence of specific growth factors, contact inhibition and anchorage dependence.


### overview: generated upcyte® cell types

We next investigated the phenotype of expanded upcyte® cells. Importantly, generated cells maintained expression of characteristic marker proteins throughout the study. For example, hepatocytes expressed cytokeratin 8/18, human serum albumin and alpha-1 antitrypsin. Accordingly, upcyte® LSECs were characterized by expression of CD31 and von-Willebrand-factor.


#### hepatocytes (HCs)




#### proximal tubular epithelial cells (PTCs)




#### liver sinusoidal endothelial cells (LSECs)




#### keratinocytes (KCs)



#### microvascular endothelial cells (mvECs)



#### skeletal muscle cells (SkMCs)




CK - cytokeratin, HSA - human serum albumin, AAT - alpha-antitrypsin, MR - mannose receptor, vWF - von-Willebrand-factor, UEA-1 - Ulex Europaeus Lectin 1, GGT1 - gamma-glutamyltransferase 1, vWF - von Willebrand Factor, eNOS - endothelial NOS or nitric oxide synthase 3 (NOS3), myosin - myosin heavy chain 2

### do you want to expand your primary cells?

#### try our upcyte® service!

Do you want to have virtually unlimited amount of your own selected donor? Do you need upcyte® cells from other cell types than what we currently offer? Do you want upcyte® cells from your diseased donor? Do you need cells from other species, e.g. monkey? Send us your primary cells - we upcyte! We will apply our upcyte® technology and produce large batches of cells from your donor of choice. The cells will be shipped back to you as cryopreserved vials or can be stored at our facility. Do not hesitate to contact us if you would like more information on how we can upcyte your favorite batch of primary cells.

### expansion strategy



#### scheme for our in house production of upcyte® cells

Starting from freshly isolated cells or a single cryopreserved vial containing  $1 \cdot 5 \cdot 10^6$  cells, cells were transduced and subsequently cultured in the presence of growth factors. Depending on the cell type used, we observed appearance of proliferating colonies at 2-3 weeks post transduction latest (e.g. hepatocytes). After the first three expansion steps, backup cells were frozen. Cells were subsequently expanded to generate a Master Cell Bank (MCB) consisting of 10 vials with  $0.5 \cdot 1 \cdot 10^7$  cells each. Single MCB vials were then used to produce a Working Cell Bank with up to 250 vials containing  $5 \cdot 1 \cdot 10^6$  cells per vial.

### examples of functional assays using upcyte® cells


#### hepatocytes: metabolism

upcyte® hepatocytes expressed metabolizing enzymes of phase I (e.g. CYP 1A2, 2B6, 2C9 and 3A4) and further exhibited phase II activities (UGT, SULT & GST). upcyte® hepatocytes further produced urea and secreted albumin (not shown). Differences in performance could be detected between cells derived from different donors.

Cells	Specific activity (pmol/min/mg protein)			
	CYP1A2	CYP2B6	CYP2C9	CYP3A4
#10-03	3.3 ± 0.4	40.3 ± 6.5	91.8 ± 5.5	21.4 ± 9.6
#151-03	0.7 ± 1.4	71.1 ± 11.3	29.1 ± 21.4	77.8 ± 22.6
#422A-03	2.3 ± 0.1	33.6 ± 11.4	4.8 ± 3.1	42.9 ± 6.3
#653-03*	17.1 ± 0.5	68.4 ± 18.4	16.2 ± 0.9	178.3 ± 17.0
HepaRG	10.0 ± 1.5	6.45 ± 0.97	4.57 ± 2.93	48.5 ± 13.9

#### LSECs: receptor-mediated endocytosis


We observed high expression of endocytic receptors, e.g. MR, FcγR and LDLR in upcyte® LSECs. Corresponding ligand uptake could be demonstrated for respective fluorophore-conjugated ligands (FITC-OVA, FITC-AGG and AF488-acLDL).



FITC-OVA = FITC ovalbumin, FITC-AGG = aggregated gamma globulin, AF488-acLDL = acetylated low density lipoprotein


#### keratinocytes: wound healing

upcyte® keratinocytes exhibit a similar migration capacity when compared to primary cells. Migration was determined by wound healing assay. Confluent monolayers of upcyte® keratinocytes were scratched using a pipette tip. Wound closure was achieved after 8 h.



#### SkMCs: atrophy & hypertrophy

Atrophy (decrease in the mass of the muscle) was induced by addition of dexamethasone & myostatin. The atrophic effect can be reversed through the inhibition with IGF-1 which also induces Hypertrophy (increase myofibril size). The hypertrophic effect can be reversed through the inhibition with LY294000



### SUMMARY & CONCLUSION

In conclusion, we developed a comprehensive platform enabling the controlled expansion of primary cells derived from various tissues for up to 40 population doublings. Importantly, upcyte® cells maintained a mature and primary-like phenotype as demonstrated by expression of marker proteins and functional assays.

We thus conclude that upcyte® expanded primary cells represent a promising model for biomedical research and drug discovery, potentially facilitating throughput and reproducibility of cell-based assays.